

## 2.1.3 Summary of methods for studying apoptosis in cell populations

Method/Roche Applied Science product	Parameter analyzed	Label	Assay principle	Advantages	Limitations	For product information, see
<b>DNA Fragmentation Assay, radioactive</b> <sup>11, 12</sup>	DNA fragmentation (LMW and HMW-DNA)	[ <sup>3</sup> H]-TdR or [ <sup>125</sup> I]-UdR, prelabel	<ul style="list-style-type: none"> <li>DNA fragments are released from the cytoplasm of apoptotic cells after lysis with non-ionic detergent.</li> <li>The LMW DNA is separated from nuclear HMW DNA by centrifugation.</li> <li>The radioactivity in the supernatant and in the pellet is determined by LSC.</li> </ul>	<ul style="list-style-type: none"> <li>Quantitative measurement over a large range (several orders of magnitude)</li> <li>Sensitive (10<sup>3</sup>-10<sup>4</sup> cells/test required)</li> <li>Suitable for analysis of cell-mediated (cytotoxicity) effects</li> </ul>	<ul style="list-style-type: none"> <li>Radioactive isotope</li> <li>Requires prelabeling and extensive washing of the target cells</li> <li>Limited to target cells proliferating <i>in vitro</i></li> <li>Increased background due to free [<sup>3</sup>H]-TdR in the cytoplasm</li> </ul>	
<b>DNA Fragmentation Assay, non-radioactive</b> <sup>13</sup> <b>Cellular DNA Fragmentation ELISA</b>	DNA fragmentation (LMW DNA)	BrdU, prelabel	<ul style="list-style-type: none"> <li>DNA fragments are released from the cytoplasm of apoptotic cells after lysis with a non-ionic detergent.</li> <li>The LMW DNA is separated from nuclear HMW DNA by centrifugation.</li> <li>The supernatant is analyzed by ELISA.</li> </ul>	<ul style="list-style-type: none"> <li>Sensitive (10<sup>3</sup>-10<sup>4</sup> cells/test required)</li> <li>Labeled cells do not have to be washed</li> <li>Optimal for microtiter plate format</li> <li>Non-radioactive</li> <li>Suitable for analysis of cell-mediated (cytotoxicity) effects</li> </ul>	<ul style="list-style-type: none"> <li>Prelabeling of the target cells required</li> <li>Can only assay target cells proliferating <i>in vitro</i></li> <li>Narrow range of quantitative measurement (only one order of magnitude)</li> </ul>	page 64 of this guide
<b>JAM Test</b> <sup>14</sup>	DNA fragmentation (HMW DNA)	[ <sup>3</sup> H]-TdR, prelabel	<ul style="list-style-type: none"> <li>Cells are harvested by vacuum aspiration onto glass fiber filters. While LMW-DNA is washed through the filters, the HMW DNA is retained on these filters.</li> <li>The radioactivity retained on the filters is measured by LSC.</li> </ul>	<ul style="list-style-type: none"> <li>Sensitive (10<sup>3</sup>-10<sup>4</sup> cells/test required)</li> <li>Only 1 washing step required for the labeled cells</li> <li>Low spontaneous release: cytotoxic events causing low cell lysis over prolonged period of time (8-24 h) can be studied</li> <li>Optimal for microtiter plate format</li> </ul>	<ul style="list-style-type: none"> <li>Radioactive isotope</li> <li>Prelabeling of the target cells required</li> <li>Limited to target cells proliferating <i>in vitro</i></li> <li>In apoptotic cells, DNA is only partially lost: viable and damaged cells are separated by only a narrow range of assay values</li> </ul>	
<b>Alkaline Elution Analysis</b> <sup>15</sup>	DNA fragmentation (LMW and HMW-DNA)	[ <sup>3</sup> H]-TdR, prelabel	<ul style="list-style-type: none"> <li>Cells are loaded onto polycarbonate filters.</li> <li>The filters are incubated with three different buffer solutions containing SDS, pH 10, SDS + Proteinase K, pH 7, or SDS, pH 12.3.</li> <li>The radioactivity in each fraction (LMW DNA) as well as the radioactivity retained on the filter (HMW DNA) is quantified by LSC.</li> </ul>	<ul style="list-style-type: none"> <li>Differential elution allows the detection of strand breaks in DNA, DNA-interstrand crosslinks and DNA-protein crosslinks</li> </ul>	<ul style="list-style-type: none"> <li>Radioactive isotope</li> <li>Prelabeling and washing of the target cells required</li> <li>Limited to target cells proliferating <i>in vitro</i></li> <li>Insensitive (10<sup>6</sup> cells/test required)</li> <li>Labor-intensive and time-consuming; only a few tests may be performed simultaneously</li> </ul>	
<b>DNA Ladder Assay</b> <sup>16</sup> (LMW and HMW DNA by size) <b>Apoptotic DNA Ladder Kit</b>	DNA fragmentation	none	<ul style="list-style-type: none"> <li>Cellular DNA is isolated by extraction and quickly purified.</li> <li>Purified total DNA (LMW and HMW DNA) is analyzed by agarose gel electrophoresis and visualized by staining with ethidium bromide.</li> </ul>	<ul style="list-style-type: none"> <li>Hallmark of apoptosis: demonstration of the mono- and oligonucleosomal DNA fragments (180 bp multimers)</li> <li>No prelabeling of the cells required: not limited to cells which proliferate <i>in vitro</i></li> <li>Non-radioactive</li> </ul>	<ul style="list-style-type: none"> <li>No quantitative measurement</li> <li>Insensitive: More than 10<sup>6</sup> cells/test required</li> <li>Labor-intensive and time-consuming; only a few tests may be performed simultaneously</li> </ul>	page 13 of this guide
<b>Nucleosome Quantification ELISA</b> <sup>13</sup> <b>Cell Death Detection ELISA</b> <sup>PLUS</sup>	DNA fragmentation (LMW DNA in association with histones)	none	<ul style="list-style-type: none"> <li>Histone complexed DNA-fragments (mono- and oligonucleosomes, LMW DNA) are released from the cytoplasm of apoptotic cells after lysis.</li> <li>The LMW DNA is separated from nuclear HMW DNA by centrifugation.</li> <li>The supernatant is analyzed by ELISA.</li> </ul>	<ul style="list-style-type: none"> <li>Sensitive (10<sup>2</sup>-10<sup>4</sup> cells/test required)</li> <li>No prelabeling of the cells required: not restricted to cells which proliferate <i>in vitro</i></li> <li>Non-radioactive</li> <li>Detection of DNA and histones in one immunoassay demonstrates mono- and oligonucleosomal DNA fragments</li> </ul>	<ul style="list-style-type: none"> <li>Samples have to be analyzed immediately because storage reduces ELISA signal</li> <li>Not recommended for tissue homogenates. Increased background could occur due to activation of nucleases during sample preparation.</li> </ul>	page 15 of this guide
<b>DNA Ladder Assay, radioactive</b> <sup>17</sup>	DNA fragmentation (LMW and HMW by size)	$\gamma$ -[ <sup>32</sup> P]-ATP, postlabel	<ul style="list-style-type: none"> <li>Cellular DNA is isolated by extraction and quickly purified.</li> <li>Purified total DNA (LMW and HMW DNA) is labeled at the 5' end with <math>\gamma</math>-[<sup>32</sup>P]-ATP by T4 Polynucleotide Kinase.</li> <li>[<sup>32</sup>P]-labeled DNA is separated by agarose gel electrophoresis and quantitated in the dried gel by a blot analyzer.</li> </ul>	<ul style="list-style-type: none"> <li>Definitive marker of apoptosis: demonstration of the mono- and oligonucleosomal DNA fragments (180 bp multimers)</li> <li>No prelabeling of the cells required: not limited to cells which proliferate <i>in vitro</i></li> <li>Highly sensitive (1000 x more sensitive than ethidium bromide): allows earlier detection of DNA fragmentation after induction of apoptosis</li> </ul>	<ul style="list-style-type: none"> <li>Labor-intensive and time-consuming; only a few tests may be performed simultaneously</li> <li>Radioactive assay (<sup>32</sup>P)</li> <li>End-labeling of purified DNA required</li> </ul>	
<b>Protease Activity Assay</b> <b>Caspase 3 Activity Assay</b>	Activation of caspases (Caspase 3)	none	<ul style="list-style-type: none"> <li>Apoptotic process including activation of the caspase cascade is induced in cells by desired method.</li> <li>Cells are lysed and cell extracts are prepared.</li> <li>Activated caspase 3 is captured out of cellular lysates by an Anti-caspase 3 antibody</li> <li>Quantification of fluorochromes cleaved from a caspase specific substrate.</li> </ul>	<ul style="list-style-type: none"> <li>Quantitative assay, cleavage of substrate is proportional to concentration of activated caspase 3 in samples</li> <li>Detection of very early stages of apoptosis</li> <li>Highly specific for caspase 3, no cross reactions with other members of the caspase family</li> </ul>	<ul style="list-style-type: none"> <li>High cell numbers needed</li> <li>Fluorescence reader, equipped with special fluorescence filters needed</li> </ul>	page 22 of this guide
<b>Protease Activity Assay</b> <b>Anti-PARP</b>	Discrete cleavage of DNA repair enzyme (PARP)	none	<ul style="list-style-type: none"> <li>Cells are treated with an apoptosis-inducing agent, which leads to induction of caspase 3 and the cleavage of Poly-ADP-Ribose-Polymerase (PARP).</li> <li>Cell extracts are prepared with SDS, fractionated by SDS-PAGE, and transferred to a PVDF membrane by western blotting.</li> <li>Blot is probed with an antibody to PARP, then with a peroxidase-labeled secondary antibody.</li> <li>Cleavage products of PARP (about 85 kD) on the membrane are revealed after an incubation with a peroxidase substrate.</li> </ul>	<ul style="list-style-type: none"> <li>Flexible, can be used with many different types of cells</li> <li>No prelabeling of cells required: not limited to cells which proliferate <i>in vitro</i></li> <li>Non-radioactive</li> <li>Marker for very early stage of apoptosis</li> </ul>	<ul style="list-style-type: none"> <li>Insensitive (requires 10<sup>5</sup>-10<sup>6</sup> cells/test)</li> <li>Labor-intensive and time-consuming; only a few tests may be performed simultaneously</li> </ul>	page 27 of this guide

▲ Table 4: Methods for studying apoptosis in cell populations.