

1.2.3 Tips for avoiding or eliminating potential TUNEL labeling artifacts

To avoid this artifact	Which may be caused by	Try the following
Nonspecific TUNEL labeling	<ul style="list-style-type: none"> ● DNA strand breaks induced by UV irradiation during tissue embedding (UV used to polymerize tissue embedding material such as methacrylate) 	<ul style="list-style-type: none"> ● Use a different embedding material, which does not require UV irradiation ● Use an alternate polymerization method
	<ul style="list-style-type: none"> ● Acid tissue fixatives (<i>e.g.</i>, mathacarn, Carnoy's fixative) cause DNA strand breaks 	<ul style="list-style-type: none"> ● Use buffered 4% paraformaldehyde as fixative
	<ul style="list-style-type: none"> ● Endogenous nuclease activity which occurs soon after tissue preparation (<i>e.g.</i>, in smooth muscle tissue slices) 	<ul style="list-style-type: none"> ● Fix tissue immediately after organ harvest ● Perfuse fixative through liver vein in intact animal
	<ul style="list-style-type: none"> ● TdT concentration too high during TUNEL labeling 	<ul style="list-style-type: none"> ● Reduce concentration of TdT by diluting it 1:2 or 1:3 with TUNEL Dilution Buffer (Cat.No. 1966 006) containing 30 mM Tris (pH 7.2) containing 140 mM sodium cacodylate and 1 mM CoCl₂
	<ul style="list-style-type: none"> ● Endogenous alkaline phosphatase activity during converter reaction 	<ul style="list-style-type: none"> ● Block endogenous AP activity by adding 1 mM levamisole to the AP substrate solution
	<ul style="list-style-type: none"> ● Endogenous peroxidase activity during converter reaction 	<ul style="list-style-type: none"> ● Before permeabilizing cells, block endogenous POD activity by immersing the slides in a solution of 0.3% H₂O₂ in methanol
	<ul style="list-style-type: none"> ● Nonspecific binding of anti-fluorescein antibody conjugate during converter reaction 	<ul style="list-style-type: none"> ● Block nonspecific sites with normal anti-sheep serum ● Block nonspecific sites with PBS containing 3% BSA (20 min) ● Use 1:2 dilution of converter solution in PBS
High background	<ul style="list-style-type: none"> ● Formalin fixation, which causes yellow staining of cells containing melanin precursors 	<ul style="list-style-type: none"> ● Use methanol fixation ● Note: <i>This fixation may lead to a reduction in TUNEL labeling sensitivity</i>
	<ul style="list-style-type: none"> ● TUNEL labeling mix too concentrated (<i>e.g.</i>, for carcinomas) 	<ul style="list-style-type: none"> ● Reduce concentration of labeling mix by diluting it 1:2 with TUNEL Dilution Buffer (Cat. No. 1966 006) containing 30 mM Tris (pH 7.2) containing 140 mM sodium cacodylate and 1 mM CoCl₂
	<ul style="list-style-type: none"> ● Endogenous alkaline phosphatase activity during converter reaction 	<ul style="list-style-type: none"> ● Block endogenous AP activity by adding 1 mM levamisole to the AP substrate solution
	<ul style="list-style-type: none"> ● Endogenous peroxidase activity during converter reaction 	<ul style="list-style-type: none"> ● Before permeabilizing cells, block endogenous POD activity by immersing the slides in a solution of 0.3% H₂O₂ in methanol
	<ul style="list-style-type: none"> ● Nonspecific binding of anti-fluorescein antibody conjugate during converter reaction 	<ul style="list-style-type: none"> ● Block nonspecific sites with normal anti-sheep serum ● Block nonspecific sites with PBS containing 3% BSA (20 min) ● Use 1:2 dilution of converter solution in PBS



To avoid this artifact	Which may be caused by	Try the following
Low TUNEL labeling (low sensitivity)	● Ethanol and methanol fixation	● Use buffered 4% paraformaldehyde as fixative
	● Extensive crosslinking during prolonged fixation reactions	● Reduce fixation time ● Use buffered 2% paraformaldehyde as fixative
	● Insufficient permeabilization of cells, so TUNEL reagents cannot reach nuclei	● Pretreat with proteinase K (concentration and time must be optimized empirically) Note: To avoid possible nuclease contamination, use only Proteinase K from Roche Applied Science, Cat. No. 03 115 836 001 ● Pretreat with 0.01 M sodium citrate for 30 min at 70°C ● Increase TUNEL incubation time
	● Restricted access of TUNEL reagents to nuclei, caused by paraffin-embedding	● After dewaxing tissue sections, treat with proteinase K (concentration, time, and temperature must be optimized empirically) Note: To avoid possible nuclease contamination, use only Proteinase K from Roche Applied Science, Cat. No 03 115 836 001 ● Immerse dewaxed tissue sections in 200 ml 0.01 M citrate buffer (pH 6.0) and treat with microwave irradiation (370 W, 5 min) Note: Conditions must be experimentally optimized for each tissue
No signal on positive control	● Inadequate DNase treatment (DNase concentration too low)	● For cryosections, apply 1 µg/ml DNase ● For paraffin-embedded tissue sections, apply 0.5 mg/ml DNase ● For many other samples, apply 1 U/ml DNase in a solution of 10 mM Tris-HCl (pH 7.4), 10 mM NaCl, 5 mM MnCl ₂ , 0.1 mM CaCl ₂ , 25 mM KCl; incubate 30 min at 37°C ● As an alternative DNase buffer, use a solution of 10 mM Tris-HCl (pH 7.5), 1 mM MgCl ₂ , 1 mg/ml BSA
Diminished TUNEL staining during DNA counterstaining	● Quenching of fluorescein signal by propidium iodide (PI)	● Use 0.5 µg/ml PI as DNA stain ● Substitute TO-PRO-3 (from Molecular Probes) in place of PI

