



▲ Figure 3: Product Selection Guide



- DNA fragmentation**
- *In Situ* Cell Death Detection Kit, POD 11 684 817 001
  - *In Situ* Cell Death Detection Kit, AP 11 684 809 001

- DNA fragmentation**
- *In Situ* Cell Death Detection Kit, TMR red 12 156 792 001
  - *In Situ* Cell Death Detection Kit, Fluorescein 11 684 795 001
  - *In Situ* Cell Death Detection Kit, POD 11 684 817 001
  - *In Situ* Cell Death Detection Kit, AP 11 684 809 001

- Membrane modification**
- Annexin-V-Biotin 10 828 690 001

- Proteolytic activity**
- M30 CytoDEATH 12 140 322 001  
12 140 349 001
  - M30 CytoDEATH, Fluorescein 12 156 857 001

- Proteolytic activity**
- M30 CytoDEATH 12 140 322 001  
12 140 349 001